

WHEN QUALITY COSTS LESS

microdot[®]

Evacuated Blood Collection Tubes




Citrate Tubes
Gel & Clot Activator Tubes
Gel & Lithium Tubes
Gel & EDTA.K2 Tubes
Glucose Tubes

Cost Effective • High Quality • Consistent draw
CE marked • Sterile • BD Compatible • Comprehensive Range


Cap Colour	IMPROVACUTER® Product	Tube Size	Type	Application
Purple	Gel & EDTA. K2 Tube	4ml (13x75mm)	EDTA	Used for viral load determination, for preparation of undiluted plasma specimens for molecular diagnostics, where undiluted plasma is required
Grey	Glucose Tube	2ml (13x75mm)	Fluoride Oxalate	Used for glucose determination
Green	Gel & Lithium Tube	6ml (13x100mm)	Lithium Heparin	Used for plasma determinations in chemistry
Yellow	Gel & Clot Activator Tube	5ml (13x100mm)	Gel	Used for serum determinations in chemistry and serology
Light Blue	Citrate Tube	2.7ml (13x100mm)	Citrate	Used for coagulation studies

Operation Flow Chart

By Multi-Sample Needle



By Blood Collection Set



Microdot Customer Care

UK: 01480 483920 / 0800 088 3920 Mon-Fri 9am-5pm ROI:1800 211 896 Mon-Fri 9am-5pm

help@microdot.biz | www.microdot.biz

microdot®

Microdot Medical Products Ltd is a wholly owned subsidiary of Cambridge Sensors Ltd.

Unit 9 Cardinal Park, Godmanchester, Huntingdon, PE29 2XG, UK

LBL-333-22A